



# Thread Lift Skin Rejuvenation

Thread Lifts by Zartaux

 Sandine  
**ZARTAU**  
DERMACEUTICAL EXPERT

# Know how:

During a thread facelift, the cosmetic surgeon inserts thin threads into the face through tiny incisions made in the targeted area. Almost immediately, treated areas will be raised slightly, eliminating sagging and giving patients a rejuvenated appearance. Thread lifting uses a series of microfilaments positioned at specific areas of the face to lift the eye brows, lift sagging skin over the cheeks, reduce nasolabial and marionette lines, define the jaw line, reducing the roundness of the face and reduce double chins. Strategic placement of the filament just underneath the skin in the direction of youthful vectors also helps stimulate collagen and firm the skin. They are 100% bio-compatible and bio-absorbable material which has long been used in medical suturing.



## Introduction:

The thread lifting treatment originated from Korea. It is a convenient, less invasive alternative to traditional facelifts, designed for people experiencing early signs of aging, such as minor facial sagging. This innovative facelift procedure can be performed in just one or two hours, often without anesthesia. If signs of aging are relatively minimal, a thread lift may be a good, or the best, way to make cosmetic enhancements while avoiding facelift scarring and other potential side effects such as downtime after the procedure or unexpected results, which sometimes occur after more intensive surgical procedures.

## Tips

Lifting action for immediate yet discreet results. This will improve over a period of 2-6 months as new collagen grows, giving further lift and tightening to the skin.


- A rejuvenating action for gradual and natural results: re-absorption of the filament component helps skin boost its natural collagen, thus further redefining the features being treated.
- The thread forms an integrated support structure for the tissue of the face due to collagen synthesis • The lifting effect is at its peak at 6 months when collagen has formed around the threads
- Lifting effect lasts from 8 months to 15 months, even though there are incidents that the results lasted up to 2 years.
- PDO threads can be safely used in combination with Botox, Fillers, Derma roller, PRP and other skin resurfacing procedures.




Mono lift or "Prime Mono" is to embed PDO into the skin. It is largely effective in facial lifting and wrinkle removal through self-stimulation of cells. There is no interruption in a daily life right after the procedure. There is no side effect or need to remove since it is biodegradable, biocompatible and dissolves in the body.



It is more advanced stage than Mono lift. It is designed as twisted pairs of the mono lift threads, resulting the acceleration of a stronger stimulation into SMAS layer and ends up forming more Collagen. Twisted PDO is thicker and stronger than the mono thread thus the effect of lifting is more powerful.



3D Gog is a barbed embedding thread. As point of contact to tow Panniculucs Carnosus is increased, Lifting effect is maximized compared to every other Lifting thread in the market.



## SCREW PDO THREAD - Collagen synthesis stimulation and wrinkles correction.

The SCREW thread is a single NORMAL thread wrapped around the needle like a spiral. The procedure using this thread looks like driving a typical screw into the tissue. This results in lifting the skin and subcutaneous tissue which causes an even greater collagen production.

### FEATURES:

- single spiral PDO thread
- skin and subcutaneous tissue lifting
- requires anaesthetic cream only, no lidocaine injection required

### TREATMENT AREAS:

- double chin correction
- eyebrows lifting
- wrinkles correction
- facial contour modelling
- body shaping
- correction of the nose shape



The SCREW PDO THREADS are available in the nine different sizes and are packed 10 pcs. of each size per box.

## Application:

### Mono / Twin

- Tightening the skin on forehead, lower lid, anterior cheek, nasolabial fold.
  1. Insert the needle.
  2. Twist at least 10 rounds to enhance more distortion for extracellular structure.
  3. Remove the needle.
    - Caution:
  4. Bleeding due to sharp needle.
  5. Compress needling sites sufficiently.

### 3D Cog

- To correct large drooping or sagging appearance.
  1. Make an inlet on skin with the needle.
  2. Insert guide wire through expected vector.
  3. Remove core wire and insert thread through external guide wire.
  4. Remove external guide wire.
  5. Pull the thread moderately to raise Cog.
  6. Cut off the stump
  - Caution: Over correction may cause dimpling or irregularity.

Apply multiple threads with minimal distortion to correct naturally.

Before

After



 Sandine  
**ZARTAUX**  
DERMACEUTICAL EXPERT